

# MEDI GUARD

MEDICAL REFRIGERATORS AND FREEZERS

## Medi Guard ULT Freezer

MG-86L1008T

# MEDI GUARD



### Touch Screen Control.

10" LCD Touch Screen controller for ease of use. Records all alarms and displays data.



### Quiet Operation.

Operates at 56.1dB(A) when compressor system is fully running. Quiet operation means less disturbance in the lab where the freezer is running.



### Security System.

Excellent alarm system equipped with a number of alarm functions - see specs for full list.



**MEDI GUARD**

Leading specialists in vaccine fridges and freezers.

[eurochill.com.au](http://eurochill.com.au)



1/13-15 Kevlar Cl, Braeside, VIC 3195



(03) 8372 0435



[sales@eurochill.com.au](mailto:sales@eurochill.com.au)

# Medi Guard ULT Freezer MG-86L1008T

PRODUCT STRUCTURE	
Volume(L)	1008
Effective Volume(L)	980.2
Upright/Chest	upright
Internal Dimensions mm(W*D*H)	1022*696*1378
External Dimensions mm(W*D*H)	1362*1025*1994
Wooden Packaging Size mm(W*D*H)	1473*1155*2176
Net weight (kg)	387
Gross weight (kg)	496
Applicable ambient temperature (°C)	16~32
Temperature setting range (°C)	-40~-86
External door type	Foam door
External material	Cold Rolled Steel Sheet
Internal material	Stainless Steel
Number of external door	2 Foam doors/handle compression
Model of temperature controller	Microprocessor
Display screen	10"HD intelligent touch screen
Drawer Material	SUS304
Shelves/Drawers	Shelves
Number of drawers	4 shelves
Size of Shelves mm	1009*649
Thermal insulation materials	VIP+ polyurethane foam
Thickness of insulation layer (mm)	door 100/side 130/back 120
Cooling mode	Direct cooling
Defrosting mode	Manual defrosting
Type of sensors	PT1000/NTC
Number of sensors	1 + 2
Number of compressors	2
Refrigerant type/grams	High Stage:R2900,150g/Low Stage:R170,125g
COMPONENTS	
Door lock	Y
Number of casters/brakes	4/2
Number of balancing feet	2
Number/inner diameter of access port	3/25mm
Data logging/Interval/Recording Time	USB/Record every 5 minutes/10years
RS485 interface	Y
Remote alarm contact	Y
High/low temperature alarm	Y
Power failure alarm	Y
Sensor fault alarm	Y
Low battery alarm	Y
Communication failure alarm	Y
Door ajar alarm	Y
PERFORMANCE	
Noise dB(A)	56.1
Load capacity of shelf (kg)	67
Temperature uniformity (°C)	1.9
Temperature fluctuation (°C)	1.6
Temperature display accuracy	0.1
Power supply (V/Hz)	220-230V, 50Hz/60Hz
Rated power (W)	850
Power consumption (kWh/24h)/(25°C ambient temperature)	10.9
Cooling rate/(min)(25°C ambient temperature)	424
Heating rate/min (25°C ambient temperature) (time to return to -50°C)	261(air)
Rated current (A)	7.44
Controller backup battery	Y/8 hours
Recovery time after 1 minute of door opening(recover to set temperature)	99

**MEDI GUARD**

Leading specialists in vaccine fridges and freezers.

[eurochill.com.au](http://eurochill.com.au)



1/13-15 Kevlar Cl, Braeside, VIC 3195



(03) 8372 0435



[sales@eurochill.com.au](mailto:sales@eurochill.com.au)